

Fusion V5 (MANE 0.95 to 1.4)

Bioinformatics Development

Sandy

2025.04.15

Transcriptome db update

- Steps (description & executed commands)
- Files summary



DB preparation steps (I)

Steps & Tools

	No.	Steps	Description	Tool
[1	download required data set from GENCODE and MANE	see sheet "Fusionv4 DB select (18,587)"	wget, rsync, zcat, samtools
2	2	generate namemap file mannually	see sheet "Fusionv4 DB select (18,587)"	awk, cat
3	3	retrieve transcript gff file	gff file preprocessing (retrieve "transcript" label from gff column 3)	zgrep, awk
		gff to bed conversion with	bed file generation (convert the information in gff to bed for transcript	convert2bed (alternative way:
4	4	"bedops_2.4.39/bin/convert2bed"	region extraction)	bedtools)
			fasta file generation (generate the fasta file for the selected regions in	
4	4	get fasta via "bedtools getfasta"	bed)	bedtools getfasta
[5 (1)	generate annotation file via RefFusion.py	<mane: 18,583=""> empty pseudo intron annotation table + pseudo N (10N) fasta generation</mane:>	RefFusion.py
į	5 (2)	generate annotation file via RefFusion.py	<gencode: 4="" curated="" id="" manually=""> empty pseudo intron annotation table + pseudo N (10N) fasta generation</gencode:>	RefFusion.py
6		ERG One manual curation step (editing annotation tables (in the empty genome + transcript loci files) + transcript fasta: keep 3 exons (in the pseudo fasta file) <- 3 output files of 5-(2)	ERG isotranscript (ENST00000442448.5) addition and adjustment	save only the first 3 exons
Ī			, , ,	,
		File concatenation with GENCODE (ERG gene adjustment)	concatenation of the curated transcript information to the exisitng	
1		and MANE output file	MANE info files obtained from previous steps (5-(1) + 6 = 7)	cat
8	3	generate annotation with GSP information	perform primer sequence alignment	blastn
8	3	generate annotation with GSP information	parse the alignment information and add it into the annotation tables	blastnparser.py

Ref repo:

https://bitbucket.org/actgenomics/reference_preparation_pipeline/src/master/



Input files/parameters (for pipeline execution)

incgctemplateFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/Blank template.fusioncalling.QC.json

Configuration

Bwa index based on fasta "refFile

```
1 => Gencode + MANE (db step 5), v0.28.0
            reffile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/y9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v0.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38 ACTFusionV4 PseudoIntron MANE-v0
v9.0 GRCh38.20231224.fasta"
          ambFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231244/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231244/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231244/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GENCODE
          annFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v0.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38 
        bwtFile = "/mnt/Bishare/sandvteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v0.95 GENCO
  9.0 GRCh38.20231224.fasta.bwt"
            pacFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v0.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v0.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38.2023124/PA043 ACTFusionV4 PseudoIntron MANE-v0.0 GRCh38.2023124/PA043 ACTFusionV4 Ps
          saFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR
v9.0 GRCh38.20231224.fasta.sa
            annoFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/PA043_ACTFusionV4_PseudoIntron_MANE-v0.95_GENCODE-r38_mPCR-v9.0_GRCh38.20231224/PA043_ACTFusionV4_PseudoIntron_MANE-v0.95_GENCODE-r38_mPCR
            gannoFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR-v9.0 GRCh38.20231224/PA043 ACTFusionV4 PseudoIntron MANE-v0.95 GENCODE-r38 mPCR
 /9.0 GRCh38.20231224.genome.5 3 UTR.loci"
           is o form fa File = "/mnt/Bishare/sandyteng/ACTFusionV4\_Torrent/data/v9.0\_input files/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa"
            isoformambFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.amb"
            isoformannFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.ann"
            isoformbwtFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.bwt"
            isoformpacFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.pac"
            isoformsaFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/y9.0 inputfiles/merged.refseg.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.sa"
            isoformmetaFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.catted.meta"
            isoformfilteringflag = "1"
            truncatedmode = "0"
           truncatedseg min aligned len = "12"
            inSpikeinFastqR1 = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/AAS502 953 IonXpress 001 1 MET 14 read-spike-in.fq.gz"
            pdbFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/protein 26db.v1.fasta"
            pdbmFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/protein_26db.v1.mapping.meta.txt"
            adapFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/torrentadapters.corrected.20230704.csv"
            acconfigFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/filter internal.QC9.0.mgsp.acr.0.5.241206.config"
            readqcconfigFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/filter fuscall 12aa Oratio read.QC.json"
```

```
boundaryqcFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/boundaryqcthresholds.v4.default.config"
    fusion container = "actgenomics/torrent fusion pipeline:v0.23.0"
process
    withLabel: process low {
        cpus = 1
        maxForks = 1
    executor = 'local'
    container = 'actgenomics/torrent fusion pipeline:v0.23.0'
    cpus = '2'
    memory = '16.G'
```

Repo: (Ref. files: ./params/fus torrent bam.14 testdata.json, ./params/fus illumina fastq.json) https://bitbucket.org/actgenomics/torrent fusion pipeline nextflow/src/master/



Modified steps from V4 → V5 transcriptome db update

DB preparation steps (V5)

• Steps & Tools

Ref. image:

• actgenomics/fusion_dev:v0.6

I	No.	Steps	Description	Tool
V5	1	download required data set from GENCODE and MANE	see sheet "Fusionv4 DB select (18,587)"=> 19,226 v1.4 MANE	wget, rsync, zcat, samtools
	2	generate namemap file manually	see sheet "Fusionv4 DB select (18,587)"=> 19,226 v1.4 MANE	awk, cat
			gff file preprocessing (retrieve "transcript" label from gff column 3)	zgrep, awk, filter_mane_gff.py
			=> Inclusion criteria: chr1-22, X, Y, protein_coding, MANE_Select	(/mnt/RD_Develop/sandyteng/
[3	retrieve transcript gff file	(summary.txt => MANE Select; manually curated)	ACTFusionV5/code/)
				convert2bed (alternative way:
			bed file generation (convert the information in gff to bed for transcript	bedtools)
4	4	gff to bed conversion with "bedops_2.4.39/bin/convert2bed"	•	(/tools/Fusion/convert2bed)
			fasta file generation (generate the fasta file for the selected regions in	
ls 4	4	get fasta via "bedtools getfasta"	bed)	bedtools getfasta
				RefFusion.py =>
				RefFusion.v2.py
			<mane: <del="">18,583> empty pseudo intron annotation table + pseudo N (10N)</mane:>	(/mnt/RD_Develop/sandyteng/
	5 (1)	generate annotation file via RefFusion.py-→ RefFusion.v2.py		ACTFusionV5/code/)
			<gencode: 4="" curated="" id="" manually=""> empty pseudo intron annotation table</gencode:>	
ŀ			+ pseudo N (10N) fasta generation	RefFusion.py
		ERG One manual curation step (editing annotation tables (in		
		the empty genome + transcript loci files) + transcript fasta:		
	_	keep 3 exons (in the pseudo fasta file) <- 3 output files of 5-		
(-	• •	, , , , , , , , , , , , , , , , , , , ,	save only the first 3 exons
		` ,	concatenation of the curated transcript information to the exisitng MANE	
	7	and MANE output file	info files obtained from previous steps (5-(1) + 6 = 7)	cat
				candidate_exons_mapping.sh
		· · · · · · · · · · · · · · · · · · ·	1039 probe regions are converted to 533 mapped (probe) exons (on	(/mnt/RD_Develop/sandyteng/
2	8-0-a	transcriptome (MANE v1.4)	pseudo-transcriptome v1.4)	ACTFusionV5/code/)
				Probe_faheader_converter.py
		modify the header of the probe fasta file (replace with		(/mnt/RD_Develop/sandyteng/
-	8-0-b	primer-like header)	Header conversion step for blastnparser.py & blastn result	ACTFusionV5/code/)
				blastn
				(/tools/Fusion/ncbi-
				blast/bin/blastn -task blastn-
1	3-0-c	generate annotation with GSP information	perform primer sequence alignment	short -dust no)
				blastnparser.py
		the state of the s		(/mnt/RD_Develop/sandyteng/
		•	parse the alignment information and add it into the annotation tables	ACTFusionV5/code/)
_	J	generate index files for bwa	generate the 5 indices required for "bwa mem"	/tools/Fusion/bwa index



Input files for probe region to MANE v1.4 loci conversion

Source files

Step 0

- Probe (region) bed file provided by Lucy
 - 1039 v5 probe regions /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/captureprobe_250401/ACTFusionv5_target-region_PartAB_individual_1039.bed
- Genome/Transcriptome DB (downloaded from websites)
 - MANE v1.4 DB (https://ftp.ncbi.nlm.nih.gov/refseq/MANE/MANE_human/release_1.4/)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/release_1.4/MANE.GRCh38.v1.4.summary.txt.gz)
 - Gencode V47 (http://ftp.ebi.ac.uk/pub/databases/gencode/Gencode_human/release_47/)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/OpenDB_GENCODE_human_r47/GRCh38.p14.genome.fa.gz)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/OpenDB_GENCODE_human_r47/gencode.v47.annotation.gff3.gz)
- Namemap
 (19,338 transcript RefSeq ID map derived from MANE v1.4 "MANE.GRCh38.v1.4.summary.txt.gz" => "MANE Select" in column 10)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.name map
- Derived files generated via steps 1-5(1) (blank loci files)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.transcript.loci (Blank Loci file)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.genome.loci (Genomic Loci file)



Source gff file

• Enter image "actgenomics/fusion_dev:v0.6" => samtools 1.9 from base image

```
./OpenDB_GENCODE_human_r47/derived/

|-- GRCh38.p14.genome.fa

|-- GRCh38.p14.genome.fa.fai

|-- gencode.v47.annotation.gff3

`-- gencode.v47.annotation.transcript.gff3
```

0 directories, 4 files

```
# commands
awk -F"\t" '{if($3=="transcript") print}' > ./OpenDB_GENCODE_human_r38/derived/gencode.v38.annotation.transcript.gff3
root@184f56381e9e:/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47# zgrep -v
"^#" ./OpenDB_GENCODE_human_r47/gencode.v47.annotation.gff3.gz | awk -F"\t" '{if($3=="transcript") print}'
> ./OpenDB_GENCODE_human_r47/derived/gencode.v47.annotation.transcript.gff3
```



Derived output summary

Derived output files

MANE v1.4 derived files

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/

- |-- MANE.GRCh38.v1.4.ensembl genomic.gff
- |-- MANE.GRCh38.v1.4.summary.namemap
- `-- MANE.GRCh38.v1.4.summary.txt

0 directories, 3 files

GENCODE v47 derived files

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/OpenDB_GENCODE_human_r47/derived/

- |-- GRCh38.p14.genome.fa
- |-- GRCh38.p14.genome.fa.fai
- -- gencode.v47.annotation.gff3
- `-- gencode.v47.annotation.transcript.gff3

0 directories, 4 files



Input files for probe region to MANE v1.4 loci conversion

- Source files
 - Probe (region) bed file provided by Lucy
 - 1039 v5 probe regions /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/captureprobe_250401/ACTFusionv5_target-region_PartAB_individual_1039.bed
 - Genome/Transcriptome DB (downloaded from websites)
 - MANE v1.4 DB (https://ftp.ncbi.nlm.nih.gov/refseq/MANE/MANE_human/release_1.4/)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/release_1.4/MANE.GRCh38.v1.4.summary.txt
 .gz)
 - Gencode V47 (http://ftp.ebi.ac.uk/pub/databases/gencode/Gencode_human/release_47/)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/OpenDB_GENCODE_human_r47/GRCh38.p14.genome.fa.gz)
 (/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/OpenDB_GENCODE_human_r47/gencode.v47.annotation.gff3.gz)

Step 2

- Namemap
 (19,338 transcript RefSeq ID map derived from MANE v1.4 "MANE.GRCh38.v1.4.summary.txt.gz" => "MANE Select" in column 10)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.name
- Derived files generated via steps 1-5(1) (blank loci files)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.transcript.loci (Blank Loci file)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.genome.loci (Genomic Loci file)



Namemap construction for 19,292 transcripts

- 66 missing IDs found (=> 66 "MANE_Plus_Clinical" IDs)
 grep -Fxv -f /mnt/RD_Develop/sandyteng/workdir/namemap_fusionv5_2_manev1.4/MANE.GRCh38.v1.4.summary.namemap.19338.namelist /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/Output_Final/MANE.GRCh38.v1.4.ensembl_genomic.transcript.list
- Generate namemap for the 66 MANE Plus Clinical transcripts

 more /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.txt | grep -wf
 /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/Output_Final/MANE.GRCh38.v1.4.summary.namemap.66.missing.namelist | awk -F"\t" '{if(\$10 == "MANE Plus Clinical")print \$8"\t"\$2"\t"\$4"\t"\$6}' >
 /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.plus.clinical.namemap
- Generate namemap for the 19338 (MANE Select) + 66 (MANE Plus Clinical) transcripts
 cat /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.namemap
 /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.plus.clinical.namemap
 /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap



GFF → bed → fasta conversion (conversion log)

/mnt/RD_Develop/sandyteng/ACTFusionV5/code/filter_mane_gff.py

```
# (1) extract gff lines with "transcript" value in the 3rd column

# (1) extract gff lines via "filter_mane_gff.py" (=> 19,292 protein coding transcripts)

python3 /mnt/RD_Develop/sandyteng/ACTFusionV5/code/filter_mane_gff.py -

i ./mane_v1.4/OpenDB_MANE_human_v1.4/release_1.4/MANE.GRCh38.v1.4.ensembl_genomic.gff.gz -

o ./mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.ensembl_genomic.transcript.gff

# (2) convert gff to bed (=> 19,292 bed lines)

/tools/Fusion/convert2bed -i gff -d < ./mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.ensembl_genomic.transcript.gff

> ./mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.ensembl_genomic.transcript.bed

# (3) extract sequence (gencode.genome.fa) via MANE v1.4 bed file (transcript.bed) (=> 19,292 fasta header lines (format: ">ENSTXXX.X"))

bedtools getfasta -name -s -fi ./gencode_v47/OpenDB_GENCODE_human_r47/derived/GRCh38.p14.genome.fa -

bed ./mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.ensembl_genomic.transcript.bed -

fo ./mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.ensembl_genomic.transcript.fasta
```



$GFF \rightarrow bed \rightarrow fasta conversion$

Convert MANE transcript gff to bed

(1) gff file filtering (includes transcript, protein coding, chr1-22, X, Y)

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane _v1.4/OpenDB_MANE_human_v1.4/derived/

- |-- MANE.GRCh38.v1.4.ensembl_genomic.gff
- |-- MANE.GRCh38.v1.4.ensembl_genomic.transcript.bed
- |-- MANE.GRCh38.v1.4.ensembl_genomic.transcript.gff
- |-- MANE.GRCh38.v1.4.summary.namemap
- `-- MANE.GRCh38.v1.4.summary.txt

0 directories, 5 files

#(2) gff to bed

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE _human_v1.4/derived/

- |-- MANE.GRCh38.v1.4.ensembl genomic.gff
- |-- MANE.GRCh38.v1.4.ensembl_genomic.transcript.bed
- |-- MANE.GRCh38.v1.4.ensembl genomic.transcript.gff
- |-- MANE.GRCh38.v1.4.summary.namemap
- `-- MANE.GRCh38.v1.4.summary.txt

0 directories, 5 files

(3) bed to fasta

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/gencode_v47/Open DB GENCODE human r47/derived/

- |-- GRCh38.p14.genome.fa
- |-- GRCh38.p14.genome.fa.fai
- -- gencode.v47.annotation.gff3
- `-- gencode.v47.annotation.transcript.gff3

0 directories, 4 files

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenD B MANE human v1.4/derived/

- |-- MANE.GRCh38.v1.4.ensembl genomic.gff
- -- MANE.GRCh38.v1.4.ensembl_genomic.transcript.bed
- -- MANE.GRCh38.v1.4.ensembl_genomic.transcript.fasta
- |-- MANE.GRCh38.v1.4.ensembl_genomic.transcript.gff
- |-- MANE.GRCh38.v1.4.summary.namemap
- `-- MANE.GRCh38.v1.4.summary.txt

0 directories, 6 files



(Rerun) Command (V5 probe → V1.4 exons) (19,292 – 66 = 19,226 MANE only transcripts)

Executed commands

generate the 19,226 id list grep -Fxv -f /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/Output_Final/MANE.GRCh38.v1.4.summary.namemap.66.missing.namelist /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/Output_Final/MANE.GRCh38.v1.4.ensembl_genomic.transcript.list > /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/Output_Final/PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0_GRCh38.20250407.transcript.MANE.only.list

```
# generate the 19,226 loci (transcript, genomic)
mkdir -p /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/Output MANE Select
python3 /mnt/RD Develop/sandyteng/ACTFusionV5/code/RefFusion.v2.py
g/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.ensembl genomic.gff
-m /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.summary.txt \
-f/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.ensembl genomic.transcript.fasta
-p /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/Output Final/PA053 ACTFusionV5 PseudoIntron MANE-v1.4 GENCODE-r47 capture-
v1.0 GRCh38.20250407.transcript.MANE.only.list
-o /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE
# create & enter container
docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion dev:v0.6
# create output folder
mkdir-p/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/
# obtain mapping exons (pseudo locations on 10*N transcriptome)
bash /mnt/RD_Develop/sandyteng/ACTFusionV5/code/candidate_exons_mapping.sh /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.genome.loc
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/Output MANE Select/20250407 MANE.transcript.loci
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/captureprobe 250401/ACTFusionv5 target-region PartAB individual 1039.bed fusionv4.MANE.v1.4.GENCODE.r47
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select//tools/Fusion
```



Get Fasta (Probe anchored exon sequences)

 Target (mapped to 1.4) bed /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/fusionv4 .MANE.v1.4.GENCODE.r47.candidate.exons.transcript.bed (bed using pseudo FASTA file generated via "RefFusion.v2.py")

```
# create output folder
mkdir -p /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseq

# create & enter container
docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion_dev:v0.6

# generate index for fasta file (for bedtools getfasta)
samtools faidx /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta

# extract mapped exons (candidate.exons.transcript.bed) sequences from gencode fasta file (gencode.genome.fa)
bedtools getfasta -name -s -fi /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta -bed
/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/fusionv4.MANE.v1.4.GENCODE.r47.candidate.exon
s.transcript.bed -fo
/mnt/RD_Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseg/MANE.GRCh38.v1.4.0407.probe.fasta
```



FASTA header conversion

- Source files
- Primer fasta file (403 primers)
 - /mnt/RD_Develop/yfhuang/ACTFusionV4DB/20240507/InhouseDB_mPCR/v9.0/20231105_mPCR9.0.fasta
- Probe fasta file (533 mapped exons)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseq/MANE_ .GRCh38.v1.4.0407.probe.fasta
- Namemap file (19,404 MANE Select + MANE Plus Clinical transcripts (MANE, v1.4))
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.se lect.and.plus.clinical.namemap
 - 19,338 MANE Select transcripts
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary.namemap
 - 66 MANE Plus Clinical transcripts
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.summary. plus.clinical.namemap



Probe faheader converter.py

usage: Probe_faheader_converter.py [-h] -f INPUTPROBEFASTA -n MANENAMEMAPTABLE -o OUTPUTPROBEFASTA

>ENST00000373237.4;PSMB2-exon-fusionv4-2()

GCCGGGCATGGTGGGCAACTGTAATCCCAGCTACTCGGGAGGCTAAGGCAAGAGAATTGCTTGAACCCAGAGGCGGAGGTTGCAGTGAACTCGAAATCGTGCCACTCCAGCCTGGG >ENST00000373237.4;PSMB2-exon-fusionv4-3()

AAACAGAGTGACTATAGTTAATAATAGCTTAACTGTATATTTTAAAATAACTAAAAGAGTATAATTGGGTT

>ENST00000373237.4;PSMB2-exon-fusionv4-5()

CGCCATTGCACTCCAGCCTGGGCAACAGAGCGAGACTCCGTCTCAAAAAA

>Probe-mane001|PSMB2|NM 002794.5|ENST00000373237.4|2|F|123

GCCGGGCATGGTGGTGGCCAACTGTAATCCCAGCTACTCGGGAGGCTAAGGCAAGAGAATTGCTTGAACCCAGAGGCGGAGGTTGCAGTGAACTCGAAATCGTGCCACTCCAGCCTGGG

>Probe-mane002|PSMB2|NM_002794.5|ENST00000373237.4|3|F|71

AAACAGAGTGACTATAGTTAATAATAGCTTAACTGTATATTTTAAAAATAACTAAAAGAGTATAATTGGGTT

>Probe-mane003|PSMB2|NM 002794.5|ENST00000373237.4|5|F|50

CGCCATTGCACTCCAGCCTGGGCAACAGAGCGAGACTCCGTCTCAAAAAA

python3 Probe_faheader_converter.py \

- -f /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseq/MANE.GRCh38.v1.4.0407.probe.fasta \
- -n/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap \
- -o MANE.GRCh38.v1.4.0407.probe.wtprimerlikeheader.fasta.gz



Step 8 (v1.4) – blastn (probe – transcriptome)

Description:

- 1. Convert Probe.fasta headers
- 2. Align probe sequence to transcriptome with pseudo-intron (10*N) via blastn

Commands

```
# probe fasta generation
python3 /mnt/RD Develop/sandyteng/ACTFusionV5/code/Probe faheader converter.py\
-f/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/probeseg/MANE.GRCh38.v1.4.0407.probe.fasta \
-n /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/mane v1.4/OpenDB MANE human v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/probeseg/MANE.GRCh38.v1.4.0407.probe.wtprimerli
keheader.fasta.gz
# unzip fasta.gz
gunzip
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/probeseg/MANE.GRCh38.v1.4.0407.probe.wtprimerli
keheader.fasta.gz
# create & enter container
docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion dev:v0.6
# folder to store blastn results
mkdir -p /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/blastn/
# add on Probe information with BLASTN result
/tools/Fusion/ncbi-blast/bin/blastn -query
/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/InhouseDB Probe/captureprobe 250407 MANE Select/probeseg/MANE.GRCh38.v1.4.0407.probe.wtprimerli
keheader.fasta -subject /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/Output MANE Select/20250407 MANE.fasta -outfmt 6 -task blastn-short >
```



Step 8 (v1.4) – blastn results

Result files

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseq/ |-- MANE.GRCh38.v1.4.0407.probe.fasta

`-- MANE.GRCh38.v1.4.0407.probe.wtprimerlikeheader.fasta

0 directories, 2 files

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/blastn/ `-- 20250407_probe.blastn

0 directories, 1 file



Step 8 (v1.4) – generate loci file with probe info

Loci file generation

create & enter container
docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion_dev:v0.6

folder to stored db output
mkdir -p /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_Loci/250407/

blastn parser

python3 /mnt/RD_Develop/sandyteng/ACTFusionV5/code/blastnparser.py \

- -if /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/blastn/20250407_probe.blastn\
- -mp /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap \
- -lf /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.transcript.loci >

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_Loci/250407/PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0 GRCh38.20250407.transcript.MANE.only.loci

ACT GENOMICS T /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_ Probe/captureprobe_250407_MANE_Select/blastn/ `-- 20250407_probe.blastn

0 directories, 1 file

`-- PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0_GRCh38.20250407.transcript.MANE.only.loci

0 directories, 1 file

727665 728813 5372 © 6520CT Lenomide CO., L5372 J. 6520 J. F. L. Proble Mango 744 served.

Step 8 (v1.4) – loci file

```
>ENST00000389048.8
                                   ENSG00000171094.18
                                                       ALK NM 004304.5
                                                                         chr2 29192774
                                                                                          29921586
                                                                                                         29
                                                                                                              728813 6520 T
                       1594 1
                                 1594 1
                                              1|1594|F|Probe-mane046
                   203890 204009 1605 1724 1
                                                    1605|1724|F|Probe-mane047
                   226573 226737 1735 1899 1
                                                    1735|1899|F|Probe-mane048
ALK (29 exons)
                   389471 389672 1910
                                      2111 1
                                                    1910|2111|F|Probe-mane049
                                                    2122|2249|F|Probe-mane050
                   537728 537855 2122 2249 1
                   593106 593237 2260
                                      2391 1
                                                    2260|2391|F|Probe-mane051
                                                                               Inspection:
                                                    2402|2533|F|Probe-mane052
                   600705 600836 2402
                                      2533 1
                                                                               Only 531 probes are annotated.
                                      2644 1
                                                    2544|2644|F|Probe-mane053
                   603183 603283 2544
                                      2824 1
                                                    2655|2824|F|Probe-mane054
                   624530 624699 2655
                    646091 646185 2835
                                       2929 1
                                                     2835 | 2929 | F | Probe-mane 055
                                                                               Missing probes:
                                                     2940|3068|F|Probe-mane056
                    646360 646488 2940
                                       3068 1
                                                                               Probe-mane081, Probe-mane501
                                                     3079|3241|F|Probe-mane057
                    670320 670482 3079
                                       3241 1
                    681757 681907 3252
                                                     3252|3402|F|Probe-mane058
                                       3402 1
                                                                               >Probe-mane501|TMPRSS2|NM 005656.4|ENST00000332149.10|2|F|71
                                                     3413|3544|F|Probe-mane059
                    687891 688022 3413
                                       3544 1
                                                                               689139 689283 3555
                                       3699 1
                                                     3555|3699|F|Probe-mane060
                                                                               GAGGGGTGAGTGAGGAGGT
                    692521 692703 3710
                                       3892 1
                                                     3710|3892|F|Probe-mane061
                    693915 694013 3903
                                       4001 1
                                                     3903 | 4001 | F | Probe-mane 062
                                                                               >Probe-mane081|ERBB4|NM 005235.3|ENST00000342788.9|7|F|142
                                       4164 1
                                                     4012 | 4164 | F | Probe-mane 063
                    694513 694665 4012
                                                                               ATATGTTATATGTTATATATATATATATGTTATATGTTATATGTTATATGTTATATATATATAT
                                                     4175 | 4279 | F | Probe-mane 064
                    696022 696126 4175
                                       4279 1
                                                                               698059 698245 4290
                                       4476 1
                                                     4290 | 4476 | F | Probe-mane 065
                                                                               TATTATATATATTATG
                                       4577 1
                    698980 699070 4487
                                                     4487 | 4577 | F | Probe-mane 066
                    699179 699243 4588
                                       4652 1
                                                     4588|4652|F|Probe-mane067
                    700752 700881 4663
                                       4792 1
                                                     4663 | 4792 | F | Probe-mane 068
                                                     4803 | 4900 | F | Probe-mane 069
                    707506 707603 4803
                                       4900 1
                    711709 711801 4911
                                       5003 1
                                                     4911|5003|F|Probe-mane070
                    714315 714416 5014 5115 1
                                                     5014|5115|F|Probe-mane071
                                                     5126|5260|F|Probe-mane072
                    723911 724045 5126 5260 1
                    724727 724817 5271 5361 1
                                                     5271|5361|F|Probe-mane073
```

The 2 missing probes

- >Probe-mane081|ERBB4|NM_005235.3|ENST00000342788.9|7|F|142
- >Probe-mane501|TMPRSS2|NM 005656.4|ENST00000332149.10|2|F|71

```
>ENST00000342788.9 ENSG00000178568.16 ERBB4 NM 005235.3 chr2 211375717
                                                                          212538802
                                                                                         28 1163086 12367 T
       354 1
               354 1
                           1|354|F|Probe-mane075
   413900 414051 365 516 1 F
                                  365|516|F|Probe-mane076
                                 527|713|F|Probe-mane077
   591187 591373 527 713 1 F
   750644 750778 724 858 1 F
                                  724|858|F|Probe-mane078
                                  869|934|F|Probe-mane079
   788099 788164 869 934 1 F
   813609 813727 945 1063 1 F
                                  945|1063|F|Probe-mane080
   816269 816410 1074 1215 0
   825155 825268 1226 1339 1
                                  1226|1339|F|Probe-mane082
                                  1350|1476|F|Probe-mane083
   826627 826753 1350 1476 1
                              F
                                                                                                  >ENST00000332149.10 ENSG00000184012.14 TMPRSS2 NM 005656.4 chr21 41464305
                                                                                                                                                                               41508158
   833412 833485 1487 1560 1 F
                                   1487 | 1560 | F | Probe-mane 084
                                                                                                  1 1 78 1 78 1 F
                                                                                                                              1|78|F|Probe-mane500
   834609 834699 1571 1661 1
                                   1571|1661|F|Probe-mane085
                                                                                                      9970 10040 89 159 0
                                   1672|1871|F|Probe-mane086
    836637 836836 1672 1871 1 F
                                                                                                                                   170|392|F|Probe-mane502
                                                                                                      13581 13803 170 392
   859619 859751 1882 2014 1
                                   1882|2014|F|Probe-mane087
                                                                                                                       489
                                                                                                                                   403 | 489 | F | Probe-mane 503
                                                                                                      18566 18652 403
                                   2025|2118|F|Probe-mane088
   865546 865639 2025 2118 1
                                                                                                      19646 19765 500 619
                                                                                                                                   500|619|F|Probe-mane504
                                   2129|2283|F|Probe-mane089
   873326 873480 2129 2283 1 F
                                                                                                                                   630|756|F|Probe-mane505
                                                                                                      27557 27683 630 756 1 F
   880975 881049 2294 2368 1
                                   2294|2368|F|Probe-mane090
                                                                                                      28877 28987 767 877 1 F
                                                                                                                                   767 | 877 | F | Probe-mane 506
   908209 908341 2379 2511 1
                                   2379|2511|F|Probe-mane091
                                                                                                                                   888|931|F|Probe-mane507
                                                                                                      31539 31582 888
                                                                                                                      931 1
                                   2522|2644|F|Probe-mane092
   914759 914881 2522 2644 1 F
                                                                                                      34663 34834 942 1113 1 F
                                                                                                                                   942 | 1113 | F | Probe-mane 508
   919528 919626 2655 2753 1
                                   2655|2753|F|Probe-mane093
                                                                                                      36178 36353 1124 1299 1 F 1124|1299|F|Probe-mane509
   976715 976900 2764 2949 1 F
                                   2764|2949|F|Probe-mane094
                                                                                                      37416 37511 1310 1405 1 F
                                                                                                                                     1310 | 1405 | F | Probe-mane 510
   1107703 1107858 2960 3115 1 F
                                    2960|3115|F|Probe-mane095
                                                                                                                                     1416 | 1558 | F | Probe-mane 511
                                                                                                      39621 39763 1416 1558 1 F
   1110320 1110395 3126 3201 1
                                    3126|3201|F|Probe-mane096
                                                                                                      40273 40425 1569 1721 1 F
                                                                                                                                     1569 | 1721 | F | Probe-mane 512
   1114502 1114648 3212 3358 1
                                    3212|3358|F|Probe-mane097
                                                                                                      42006 43854 1732 3580 1 F 1732 3580 F Probe-mane 513
                                    3369|3466|F|Probe-mane098
   1116699 1116796 3369 3466 1
   1118192 1118362 3477 3647 1 F
                                    3477|3647|F|Probe-mane099
                                     3658|3705|F|Probe-mane100
   1150811 1150858 3658 3705 1
   1151653 1151950 3716 4013 1
                               F
                                    3716|4013|F|Probe-mane101
                                     4024|12367|F|Probe-mane102
   1154743 1163086 4024 12367 1 F
```



Try different blastn command (to rescue the 2 missing probes)

Test commands => the 2 missing probes are rescued!

docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion_dev:v0.6

add on Probe information with BLASTN result
/tools/Fusion/ncbi-blast/bin/blastn -query
/mnt/RD_Develop/sandyteng/ACTEusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseg/MANE_GRCh38.v1.4.04

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/probeseq/MANE.GRCh38.v1.4.0407.probe.wtprimerlikehe ader.fasta -subject /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta -outfmt 6 -task blastn-short -dust no /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/blastn/20250407_probe.default.blastn

blastn parser (loci annotation)
python3 /mnt/RD_Develop/sandyteng/ACTFusionV5/code/blastnparser.py \
-if /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/blastn/20250407_probe.default.blastn \
-mp /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/mane_v1.4/OpenDB_MANE_human_v1.4/derived/MANE.GRCh38.v1.4.select.and.plus.clinical.namemap \
-lf /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.transcript.loci >
/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_Loci/250407/PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0_GRCh38.20250407.transcript.MANE.only.blastn.default.loci



create & enter container

BWA indices generation

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select

- |--/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta
- |--/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta.fai
- |--/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.genome.loci
- $|--/mnt/RD_Develop/s and yteng/ACTF usion V5/db_fusion v5/Output_MANE_Select/20250407_MANE. is otranscript_retrieval.log$
- `--/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.transcript.loci

0 directories, 5 files

- Source FASTA (19,226 MANE Select v1.4 transcripts)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta

create & enter container docker run --rm -v /mnt/:/mnt/ -it actgenomics/fusion dev:v0.6

generate bwa indices for fusion calling module ("ACTGfuscall.py")
/tools/Fusion/bwa index /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta

```
oot@d550999aal9b:/tools# /tools/Fusion/bwa index /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output MANE Select/20250407 MANE.fasta
[bwa index] Pack FASTA... 0.84 sec
[bwa index] Construct BWT for the packed sequence...
[BWTIncCreate] textLength=145800268, availableWord=22258764
[BWTIncConstructFromPacked] 10 iterations done. 36716444 characters processed.
[BWTIncConstructFromPacked] 20 iterations done. 67829580 characters processed.
[BWTIncConstructFromPacked] 30 iterations done. 95479004 characters processed.
[BWTIncConstructFromPacked] 40 iterations done. 120049852 characters processed.
[BWTIncConstructFromPacked] 50 iterations done. 141884412 characters processed.
[bwt gen] Finished constructing BWT in 52 iterations.
[bwa index] 35.79 seconds elapse.
[bwa_index] Update BWT... 0.86 sec
[bwa index] Pack forward-only FASTA... 0.57 sec
[bwa index] Construct SA from BWT and Occ... 20.03 sec
[main] Version: 0.7.17-r1188
[main] CMD: /tools/Fusion/bwa index /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta
 main] Real time: 61.155 sec; CPU: 58.087 sec
```

output directory

/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/

|-- 20250407 MANE.fasta

|-- 20250407_MANE.fasta.amb

|-- 20250407 MANE.fasta.ann

|-- 20250407 MANE.fasta.bwt

|-- 20250407_MANE.fasta.fai

|-- 20250407 MANE.fasta.pac

-- 20250407_MANE.fasta.sa

|-- 20250407 MANE.genome.loci

|-- 20250407 MANE.isotranscript retrieval.log

`-- 20250407 MANE.transcript.loci

0 directories, 10 files



File summary

- Steps
- Generated files
- Updated configuration file



DB preparation steps (V5) (clean)

Steps & Tools

Ref. image:

• actgenomics/fusion_dev:v0.6

	No.	Steps	Description	Tool
	1	download required data set from GENCODE and MANE	see sheet "Fusionv4 DB select (18,587)"=> 19,226 v1.4 MANE	wget, rsync, zcat, samtools
S	2	generate namemap file manually	see sheet "Fusionv4 DB select (18,587)"=> 19,226 v1.4 MANE	awk, cat
3			gff file preprocessing (retrieve "transcript" label from gff column 3)	zgrep, awk, filter_mane_gff.py
			=> Inclusion criteria: chr1-22, X, Y, protein_coding, MANE_Select	(/mnt/RD_Develop/sandyteng/
	3	retrieve transcript gff file	(summary.txt => MANE Select; manually curated)	ACTFusionV5/code/)
				convert2bed (alternative way:
			bed file generation (convert the information in gff to bed for transcript	bedtools)
	4	gff to bed conversion with "bedops_2.4.39/bin/convert2bed"	region extraction)	(/tools/Fusion/convert2bed)
			fasta file generation (generate the fasta file for the selected regions in	
	4	get fasta via "bedtools getfasta"	bed)	bedtools getfasta
				RefFusion.py =>
				RefFusion.v2.py
			<mane: 18,583=""> empty pseudo intron annotation table + pseudo N (10N)</mane:>	(/mnt/RD_Develop/sandyteng/
	5 (1)	generate annotation file via RefFusion.py → RefFusion.v2.py	fasta generation => 19,226 v1.4 MANE	ACTFusionV5/code/)
				candidate_exons_mapping.sh
		convert v5 probe regions to the regions on pseudo	1039 probe regions are converted to 533 mapped (probe) exons (on	(/mnt/RD_Develop/sandyteng/
	8-0-a	transcriptome (MANE v1.4)	pseudo-transcriptome v1.4)	ACTFusionV5/code/)
				Probe_faheader_converter.py
		modify the header of the probe fasta file (replace with		(/mnt/RD_Develop/sandyteng/
	8-0-b	primer-like header)	Header conversion step for blastnparser.py & blastn result	ACTFusionV5/code/)
				blastn
				(/tools/Fusion/ncbi-
				blast/bin/blastn -task blastn-
	8-0-c	generate annotation with GSP information	perform primer sequence alignment	short -dust no)
				blastnparser.py
				(/mnt/RD_Develop/sandyteng/
	8-0-c	generate annotation with GSP information	parse the alignment information and add it into the annotation tables	ACTFusionV5/code/)
	9	generate index files for bwa	generate the 5 indices required for "bwa mem"	/tools/Fusion/bwa index
	5	generate index files for bwa	generate the 3 maices required for bwa mem	/ tools/ Fusion/ bwa muex



Generated files summary (clean)

Probe meta file (533 mapped (probe) exons)

Step 8-0-a

- /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/InhouseDB_Probe/captureprobe_250407_MANE_Select/fusionv4.MANE.v1.4.GENCODE.r47.can didate.exons.transcript.bed
- Transcriptome (with 10*N)
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta (FASTA file for the 19,226 selected transcripts)

Step 5 (1)

- /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_MANE_Select/20250407_MANE.fasta.
bwa index suffix> Step 9 (suffices include "sa", "pac", "ann", "amb" and "bwt")
- Transcript Loci files
 - /mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/Output MANE Select/20250407 MANE.transcript.loci (Blank Loci file) Step 5 (1)

Step 8-0-c

- /mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/Output_Loci/250407/PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0 GRCh38.20250407.transcript.MANE.only.blastn.default.loci (Loci file with 533 probes (mapped exons in v1.4) annotated)
- Genomic Loci file



Updated pipeline input files (colored)

/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/

0 directories, 10 files



Updated configuration file

Configuration file (for db update test)

```
refFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db_fusionv5/v1.4_inputfiles/20250407_MANE.fasta"
   ambFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/20250407 MANE.fasta.amb'
   annFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/20250407 MANE.fasta.ann'
   bwtFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/20250407 MANE.fasta.bwt'
   pacFile = "/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/v1.4_inputfiles/20250407_MANE.fasta.pac
   saFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/20250407 MANE.fasta.sa"
    annoFile = "/mnt/RD_Develop/sandyteng/ACTFusionV5/db_fusionv5/v1.4_inputfiles/PA053_ACTFusionV5_PseudoIntron_MANE-v1.4_GENCODE-r47_capture-v1.0_GRCh38.20250407.transcript.MANE.only.blastn.default.loci
    gannoFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/20250407 MANE.genome.loci"
    isoformfaFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseg.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa"
    isoformambFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseg.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.amb"
    isoformannFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.ann"
    isoformbwtFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.bwt"
    isoformpacFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.pac"
    isoformsaFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseg.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.sa"
    isoformmetaFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/merged.refseq.20240807.gencode.v38.transcripts.protein.coding.22.isoforms.fa.catted.meta"
    isoformfilteringflag = "1"
   truncatedmode = "0"
   truncatedseq min aligned len = "12"
    inSpikeinFastqR1 = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/AAS502 953 IonXpress 001 1 MET 14 read-spike-in.fq.gz"
    pdbFile = "/mnt/Bishare/sandyteng/ACTFusionV4_Torrent/data/v9.0_inputfiles/protein_26db.v1.fasta"
    pdbmFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/protein 26db.v1.mapping.meta.txt"
    adapFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/torrentadapters.corrected.20230704.csv"
    qcconfigFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/filter internal.QC9.0.mgsp.qcr.0.5.241206.config"
   readqcconfigFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/filter fuscall 12aa Oratio read.QC.ison"
    primerlabelFile = "/mnt/RD Develop/sandyteng/ACTFusionV5/db fusionv5/v1.4 inputfiles/fusionv4.MANE.v1.4.GENCODE.r47.candidate.exons.transcript.bed
    incgctemplateFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/Blank template.fusioncalling.QC.json"
                       boundaryqcFile = "/mnt/Bishare/sandyteng/ACTFusionV4 Torrent/data/v9.0 inputfiles/boundaryqcthresholds.v4.default.config"
    fusion container = "actgenomics/torrent fusion pipeline:v0.23.0"
process {
    withLabel: process low
        cpus = 1
        maxForks = 1
   executor = 'local'
   container = 'actgenomics/torrent fusion pipeline:v0.23.0'
    memory = '16.G'
```



v1.4 annotation evaluation

- Result summary
- Manual inspection



Date: 2025-4-14 (done)

Fusionv4 test v0.1 (config test: MANE v0.95 \rightarrow v1.4)

300x IVTRNA

```
sandyteng@tp-rd177:/mmt/RD Develop/sandyteng/ACTFusionV5/nextflow outdir/20250414 ACTFusion v4 81-IVTRNA-240906-300x v0.1 fusionv4 MANE v1.4 config test$ 1s -al
total 8
drwxrwxr-x 2 sandyteng sandyteng 4096 Apr 14 15:12 .
drwxrwxr-x 9 sandyteng sandyteng 4096 Apr 14 15:12 ..
sandyteng@tp-rdl77:/mnt/RD Develop/sandyteng/ACTFusionV5/nextflow outdir/20250414 ACTFusion v4 81-IVTRNA-240906-300x v0.1 fusionv4 MANE v1.4 config test$ nextflow run /mnt/RD Develop
peline nextflow/Workflow docker illumina.nf \
-params-file /mnt/RD Develop/sandyteng/ACTFusionV5/nextflow/repo code vl.4 dbtest 0414.2025/params/ACTFusion v4 81-IVTRNA-240906-300x v0.1 fusionv4 MANE vl.4 test.json \
-c /mnt/RD_Develop/sandyteng/ACTFusionV5/nextflow/repo_code_v1.4_dbtest_0414.2025/dockerconfigs/fusion_multi_localdocker.v9.20241125.v0.23.0_v1.4.MANE.transcriptome.config
 extflow 24.10.5 is available - Please consider updating your version to it
NEXTFLOW
                  ~ version 24.04.4
Launching `/mnt/RD Develop/sandyteng/ACTFusionV5/repo dev/hybridcapture fusion pipeline nextflow/Workflow docker illumina.nf` [modest venter] DSL2 - revision: ecbd2c7c06
         e] process > SampleInfo (ACTFusion v4 81-IVTRNA-240906-300x R1.fastq.gz ACTFusion v4 81-IVTRNA-240906-300x R2.fastq.gz)
                                                                                                                                                 [100%] 1 of 1 v
   4be765] process > Fusion:mergefastq (ACTFusion v4 81-IVTRNA-240906-300x ACTFusion v4 81-IVTRNA-240906-300x R1.fastq.gz)
                                                                                                                                                 [100%] 1 of 1 v
         0] process > Fusion:trimadap (ACTFusion v4 81-IVTRNA-240906-300x ACTFusion v4 81-IVTRNA-240906-300x.fastp.zcatted.merged.fastq.gz)
                                                                                                                                                 [100%] 1 of 1 v
   /15f2ce] process > Fusion:fastp (ACTFusion v4 81-IVTRNA-240906-300x ACTFusion v4 81-IVTRNA-240906-300x.fastp.zcatted.merged .trimadap.fq.gz) [100%] 1 of 1 🔻
         process > Fusion:bwa isoform (ACTFusion v4 81-IVTRNA-240906-300x ACTFusion v4 81-IVTRNA-240906-300x.clean.fg.gz)
                                                                                                                                                 [100%] 1 of 1 •
         s] process > Fusion:bwa se (ACTFusion v4 81-IVTRNA-240906-300x ACTFusion v4 81-IVTRNA-240906-300x.isofilered.fq.gz)
                                                                                                                                                 [100%] 1 of 1 v
  /fc33la] process > Fusion:bwa fusioncalling (ACTFusion v4 81-IVTRNA-240906-300x)
                                                                                                                                                 [100%] 1 of 1 v
   [680b36] process > Fusion:fuscall2QC illumina (ACTFusion v4 81-IVTRNA-240906-300x)
                                                                                                                                                 [100%] 1 of 1 v
 ompleted at: 14-Apr-2025 15:20:22
Duration : 5m 26s
CPU hours : 0.2
Succeeded : 8
```

nextflow run /mnt/RD_Develop/sandyteng/ACTFusionV5/repo_dev/hybridcapture_fusion_pipeline_nextflow/Workflow_docker_illumina.nf \
-params-file /mnt/RD_Develop/sandyteng/ACTFusionV5/nextflow/repo_code_v1.4_dbtest_0414.2025/params/ACTFusion_v4_81-IVTRNA-240906300x v0.1 fusionv4 MANE v1.4 test.json \

-c/mnt/RD_Develop/sandyteng/ACTFusionV5/nextflow/repo_code_v1.4_dbtest_0414.2025/dockerconfigs/fusion_multi_localdocker.v9.20241125.v0.23.0_v1.4.MANE.transcriptome.config

Result summary (v1.4 annotation evaluation)

- Test sample: ACTFusion v4 81-IVTRNA-240906-300x
 - 300 replicates of the 81 IVT sequences
 - Summarized result
 - /mnt/RD_Develop/sandyteng/ACTFusionV5/nextflow/repo_code_v1.4_dbtest_0414.2025/summary/
 - Raw result folder
 - /mnt/RD_Develop/sandyteng/ACTFusionV4/nextflow_outdir/20241121_v0.20.0_image_v0.25.0_pipeline_test/ACTFusion_v4_81-IVTRNA-240906_check/ACTFusion_v4_81-IVTRNA-240906/
- Forward primers successfully annotated (for the 23 identified fusions & 3 exon-skipping cases)
- 81-read inspection
 - 17 missing reads
 - Not recognized by "fuscall.py" program, but primary alignments of the 17 reads are found within the aligned.bwa.bam (unaligned reads => SAM bitwise FLAG == 4)
 - 64 reads processed by "ACTGfuscall.py"
 - 36 WILDTYPE (inclusive of AR-V7) (22 wild type reads with no primer detected)
 - 23 FUSION
 - 3 EXONSKIPPING (MET, EGFR (2), BRAF)
 - 2 AMBIGUITY (FusionRef_659, FusionRef_689)
- Recall: 23/81 ~ 28.4%
 - (81-23)/81 ~ 2/3 of the IVT variants designed using v0.95 transcripts cannot be detected when migrating the transcriptome to v1.4 (which is reasonable since v1.4 and v0.95 included different preferred transcripts)



Result (list of missing fusion reads)

17 missing reads

IVT-RNA ID	Boundary	Туре	Group	Report status (v0.24.0)
FusionRef_077	CD74:6-ROS1:33	FUSION	1	+
FusionRef_640	AKT3:12-HEATR1:30	FUSION	1	+
FusionRef_680	SRGAP3:7-RAF1:8	FUSION	1	+
FusionRef_646	ERG:6-APP:14	FUSION	2	+
FusionRef_701	EZR:10-ROS1:35	FUSION	2	+
FusionRef_707	KIAA1549:15-BRAF:9	FUSION	2	+
FusionRef_706	KIAA1549:16-BRAF:9	FUSION	3	+
FusionRef_652	NRG1:7-PSAP:11	FUSION	4	+
FusionRef_699	CD74:6-ROS1:35	FUSION	4	+
FusionRef_704	TMPRSS2:2-ERG:2	FUSION	5	+
FusionRef_669	EZR:12-ERBB4:18	FUSION	6	+
FusionRef_709	FGFR2:17-SHTN1:7	FUSION	6	+
FusionRef_662	RSPO2:4-EIF3E:3	FUSION	7	+
FusionRef_683	EIF3E:2-RSPO2:4	FUSION	7	+
FusionRef_684	FKBP5:6-TMPRSS2:7	FUSION	7	+
FusionRef_677	SIL1:3-NRG2:6	FUSION	8	+



Result (23 FUSION IVT reads)

23 Fusions (20 of which are equivalent to the boundaries reported via MANE v0.95)

			-	_							2.5	
Fusion Type	T_Type	G_Type	5-gsp_name	5-gsp_coverage	5-1 ranscript	5-orientation	5-Exon	3-gsp_name	3-gsp_coverage 3-Transcript	3-orientation	3-Exon	Boundary Check
FUSION	ENST00000292431.5:5-ENST00000277120.8:11	NACC1:5-NTRK2:11	-	-	-	-	-	Probe-mane376	0.574257ENST00000277120.8	F	11	FALSE
FUSION	ENST00000263966.8:4-ENST00000263967.4:15	USP13:4-PIK3CA:15	-	-	-	-	-	Probe-mane118	0.495327ENST00000263967.4	F	15	TRUE
FUSION	ENST00000240851.9:4-ENST00000524377.7:9	TFG:4-NTRK1:9	-	-	-	-	-	Probe-mane038	0.910714 ENST00000524377.7	F	10	TRUE
FUSION	ENST00000382051.8:4-ENST00000397752.8:15	SLC34A2:4-MET:15	-	-	-	-	-	Probe-mane309	0.251082 ENST00000397752.8	F	15	TRUE
FUSION	ENST00000269571.10:23-ENST00000447079.6:6	ERBB2:23-CDK12:6	Probe-mane479	0.408163265	ENST00000269571.10	F	23	} -	-	-	-	TRUE
FUSION	ENST00000206249.8:3-ENST00000371998.8:15	ESR1:3-NCOA3:15	Probe-mane243	0.478632479	ENST00000206249.8	F	3	 -	-	-	-	TRUE
FUSION	ENST00000288135.6:20-ENST00000257290.10:6	KIT:20-PDGFRA:6	Probe-mane164	0.537735849	ENST00000288135.6	F	20)-	-	-	-	TRUE
FUSION	ENST00000397752.8:20-ENST00000358204.9:3	MET:20-TES:3	Probe-mane314	0.416058394	ENST00000397752.8	F	20)-	-	-	-	TRUE
FUSION	ENST00000524377.7:16-ENST00000643579.2:7	NTRK1:16-TPM4:7	Probe-mane044	0.352201258	ENST00000524377.7	F	16	j-	-	-	-	FALSE
FUSION	ENST00000277120.8:18-ENST00000314328.14:8	NTRK2:18-SPECC1L:8	Probe-mane383	0.389937107	ENST00000277120.8	F	18	} -	-	-	-	TRUE
FUSION	ENST00000338663.12:2-ENST00000405005.8:6	SLC3A2:2-NRG1:6	-	-	-	-	-	Probe-mane345	0.384615 ENST00000405005.8	F	6	TRUE
FUSION	ENST00000643579.2:7-ENST00000524377.7:10	TPM4:7-NTRK1:10	-	-	-	-	-	Probe-mane038	0.946429 ENST00000524377.7	F	10	FALSE
FUSION	ENST00000508803.6:5-ENST00000440486.8:10	NSD2:5-FGFR3:10	-	-	-	-	-	Probe-mane175	0.376712 ENST00000440486.8	F	10	TRUE
FUSION	ENST00000440486.8:17-ENST00000313288.9:8	FGFR3:17-TACC3:8	Probe-mane182	0.481132075	ENST00000440486.8	F	17	'-		-	-	TRUE
FUSION	ENST00000545606.6:8-ENST00000275493.7:15	CAND1:8-EGFR:15	-	-	-	-	-	Probe-mane263	0.35443 ENST00000275493.7	F	15	TRUE
FUSION	ENST00000323571.9:5-ENST00000269571.10:13	WIPF2:5-ERBB2:13	-	-	-	-	-	Probe-mane469	0.571429ENST00000269571.1	0 F	13	TRUE
FUSION	ENST00000355710.8:19-ENST00000163416.7:4	RET:19-GOLGA5:4	Probe-mane421	0.398648649	ENST00000355710.8	F	19)-	-	-	-	TRUE
FUSION	ENST00000440486.8:17-ENST00000313288.9:11	FGFR3:17-TACC3:11	Probe-mane182	0.481132075	ENST00000440486.8	F	17	'-	-	-	-	TRUE
FUSION	ENST00000361752.8:6-ENST00000277120.8:14	QKI:6-NTRK2:14	-	-	-	-	-	Probe-mane379	0.301587ENST00000277120.8	F	14	TRUE
FUSION	ENST00000440486.8:17-ENST00000313288.9:10	FGFR3:17-TACC3:10	Probe-mane182	0.481132075	ENST00000440486.8	F	17	'-	-	-	-	TRUE
FUSION	ENST00000396373.9:4-ENST00000277120.8:14	ETV6:4-NTRK2:14	-	-	-	-	-	Probe-mane379	0.301587ENST00000277120.8	F	14	TRUE
FUSION	ENST00000431975.7:5-ENST00000206249.8:5	AKAP7:5-ESR1:5	-	-	-	-	-	Probe-mane245	0.42446 ENST00000206249.8	F	5	TRUE
FUSION	ENST00000264312.12:8-ENST00000288135.6:8	OCIAD1:8-KIT:8	-	-	-	-	-	Probe-mane152	0.469565 ENST00000288135.6	F	8	TRUE



Result (FUSION IVT read)

• NTRK1:16-TPM3:8 (MANE v0.95) \rightarrow NTRK1:16-TPM4:7 (MANE v1.4)

MANE v0.95

FusionRef_654 FUSION ENST00000524377.7:16-ENST00000651641.1:8 NTRK1:16-TPM3:8 ENST00000524377.7:16-ENST00000651641.1:8,9 NTRK1:16-TPM3:8,9 1 56
TCTGGGAGATCTTCACCTACGGCAAGCAGCCCTGGTACCAGCTCTCCAACACGGAGGCAGAGACCCGTGCTGAGGAGATCGGTAGCCAAGCTGGAAAAGACAATTGATGACCTGGAAGATGAGC
TCTATGCCCAGAAACTGAAGTACAAGGCCATTAGCGAGGAGCTGGAC FT-mPCR0022 1 ENST00000524377.7 F 16 - - - - - - - -

MANE v1.4

FusionRef_654 FUSION ENST00000524377.7:16-ENST00000643579.2:7 NTRK1:16-TPM4:7 ENST00000524377.7:16-ENST00000643579.2:7 NTRK1:16-TPM4:7 1 57
TCTGGGAGATCTTCACCTACGGCAAGCAGCCCTGGTACCAGCTCTCCAACACGGAGGCAGAGCCCGTGCTGAGTTTGCTGAGAGATCGGTAGCCAAGCTGGAAAAGACAATTGATGACCTGGAAGATGAGC
TCTATGCCCAGAAACTGAAGTACAAGGCCATTAGCGAGGAGCTGGAC Probe-mane044 0.3522012578616352 ENST00000524377.7 F 16 - - - - - -

Genomic coordinates

MANE v0.95

NTRK1(NM_002529.4), chr1:156880157; TPM3(NM_152263.4), chr1:154170469

MANE v1.4

NTRK1(NM_002529.4), chr1:156880157; TPM4(NM_003290.3), chr19:16093684



Result (3 EXONSKIPPING IVT reads)

- 3 (/ 5 EXONSKIPPING + ARV7) targets are reported
 - BRAF:19-BRAF:11 are not reported

IVT-RNA ID	Boundary	Туре	Group	Report status (v0.24.0)	Fusion Type	G_Type	5-gsp_name	5-gsp_coverage	5-Transcript	5-orientation	5-Exon	3-gsp_name	3-gsp_coverage	3-Transcript	3-orientation	3-Exon
ARV7	AR:2,3,4	WILDTYPE	9	+	WILDTYPE	AR:2	Probe-mane527	0.151315789	ENST00000374690.9	F	2	Probe-mane528	1	ENST00000374690.9	F	3
FusionRef_685	MET:13-MET:15	EXONSKIPPING	9	+	EXONSKIPPING	MET:13-MET:15	Probe-mane307	0.394904459	ENST00000397752.8	F	13	Probe-mane309	0.251082251	ENST00000397752.8	F	15
FusionRef_686	EGFR:1-EGFR:8	EXONSKIPPING	9	+	EXONSKIPPING	EGFR:1-EGFR:8	Probe-mane249	0.166189112	ENST00000275493.7	F	1	Probe-mane256	0.452991453	ENST00000275493.7	F	8
FusionRef_687	BRAF:19-BRAF:11	EXONSKIPPING	9	+												
FusionRef_688	EGFR:24-EGFR:18	EXONSKIPPING	9	+	EXONSKIPPING	EGFR:24-EGFR:18	Probe-mane272	0.642857143	ENST00000275493.7	F	24	Probe-mane266	0.504065041	ENST00000275493.7	F	18



Result (WILDTYPE inspection, FusionRef_665)

Fusion V4 (MANE v0.95) & Arriba calling results

NM IDs reported by Arriba for the same boundary

MAGI3(NM_001142782	AKT3(NM 181690)
100 (0.00)	, , , , , , , , , , , , , , , , , , , ,

Boundary;(5' gene coordinate,3' gene coordinate)	Boundary	Туре	Group	IVT-RNA ID	Report status (v0.24.0)	5' NM ID	3' NM ID	5' gene coordinate	3' gene coordinate
MAGI3-AKT3;(chr1:113622994,chr1:243695716)	MAGI3:9-AKT3:3	FUSION	6	FusionRef_665	+	MAGI3(NM_001142782.2)	AKT3(NM_005465.7)	chr1:113622994	chr1:243695716



Result (WILDTYPE inspection, FusionRef 665)

 V5 (MANE v1.4) Result / blastn result

Note:

- 1. Different sets of preferred transcriptome => different alignment result
- 2. If we includes multiple transcripts for the same gene, it may result in ambiguity in terms of the alignment result and the corresponding NM ID(s) to report.

Calling result

FusionRef 665 WILDTYPE

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value ▼	Per. Ident	Acc. Len	Accession
✓	$\underline{\text{Homo sapiens membrane associated } \underline{\text{guanylate kinase}}, \underline{\text{WW and PDZ domain containing 3}} (\underline{\text{MAGI3}}), \underline{\text{tran}}}$	Homo sapiens	257	257	67%	4e-64	100.00%	7559	NM_152900.3
	Homo sapiens membrane-associated guanylate kinase-related MAGI-3 (MAGI-3) mRNA, complete cds	Homo sapiens	257	257	67%	4e-64	100.00%	3504	AF213259.1
✓	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	6965	XM_054335878.1
✓	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	6965	XM_005270737.4
	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	3473	XM_054335879.1
	Homo sapiens membrane-associated guanylate kinase MAGI3 mRNA, complete cds	Homo sapiens	257	257	67%	4e-64	100.00%	4288	AF257238.1
	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	6993	XM_054335877.1
	Human DNA sequence from clone RP11-473L1 on chromosome 1, complete sequence	Homo sapiens	257	257	67%	4e-64	100.00%	71464	AL390759.10
	Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3, mRNA (cDNA	Homo sapiens	257	257	67%	4e-64	100.00%	3580	BC130409.1
	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	6993	XM_017000974.2
	Homo sapiens cDNA clone IMAGE:9052664, containing frame-shift errors	Homo sapiens	257	257	67%	4e-64	100.00%	3653	BC144144.1
✓	PREDICTED: Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3	Homo sapiens	257	257	67%	4e-64	100.00%	3473	XM_047417371.1
~	Homo sapiens membrane associated guanylate kinase, WW and PDZ domain containing 3 (MAGI3), tran	Homo sapiens	257	257	67%	4e-64	100.00%	6874	NM_001142782.2
0000	00007F4C 14.0 NAACI2.0 1								

ENST00000307546.14:9 MAGI3:9 1

GGGTTGTTCGAGCATCACTGAAAAAAAGCACAATGGGATTTGGTTTTACTATTATTGGTGGAGATAGACCTGATGAGTTCCTACAAGTGAAAAAATGTGCTGAAAGATGGTCCCGCAGCTCAGGATG

MANE (RefSeq – Ensembl id map)

ENST00000307546.14 ENSG00000081026.19

MAGI3 NM 001142782.2 # the same id as the one in fusion v4 (MANE v0.95)

AKT3 NM 005465.7 # the same id as the one in fusion v4 (MANE v0.95) ENST00000673466.1 ENSG00000117020.19

Boundary;(5' gene coordinate,3' gene coordinate)	Boundary	Туре	Group	IVT-RNA ID	Report status (v0.24.0)	5' NM ID	3' NM ID	5' gene coordinate	3' gene coordinate
MAGI3-AKT3;(chr1:113622994,chr1:243695716)	MAGI3:9-AKT3:3	FUSION	6	FusionRef_665	+	MAGI3(NM_001142782.2)	AKT3(NM_005465.7)	chr1:113622994	chr1:243695716

